# DRD project

Lui Maria

March 14, 2020

## 1 Steps of the analysis pipeline

1. Load raw data with minfi and create an object called RGset storing the RGChannelSet object

```
getwd()
 # Load idat files and the experiment samplesheet using minfi
3 library (minfi)
4 vignette("minfi")
 # Set the directory in which the raw data are stored
6 baseDir <- ("/home/maria/Desktop/DRD2/Input_data")
 targets <- read.metharray.sheet(baseDir)</pre>
 targets
10 # Create an object of class RGChannelSet
RGset <- read.metharray.exp(targets = targets)
12 #These classes represents raw (unprocessed) data from a two
     color micro array; specifically an Illumina methylation array
_{14} # Now I change the directory in which I want to save the new
     files
15 setwd("/home/maria/Desktop/DRD2/newfiles")
save (RGset, file="RGset.RData")
 RGset
18
19 ? RGChannelSet
```

2. Create the dataframes Red and Green to store the red and green fluorescence respectively .

```
Red <- data.frame(getRed(RGset))
#dim(Red)
Green <- data.frame(getGreen(RGset))
#dim(Green)</pre>
```

SAMPLE	RED	GREEN	TYPE
X5775278008_R01C01	3445	858	II
X5775278008_R02C01	5408	1223	II
X5775278008_R04C01	863	1204	II
X5775278008_R05C01	1053	1044	II
X5775278035_R01C01	2114	505	II
X5775278035_R02C01	468	579	II
X5775278035_R04C01	6395	1234	II
X5775278035_R05C01	1292	1120	II
X9344737127_R01C02	330	354	II
X9344737127_R03C02	282	493	II
X9344737127_R06C01	330	468	II
X9344737127_R02C02	379	421	II
X9376538140_R02C01	399	390	II
X9376538140_R05C01	333	577	II
X9376538140_R01C01	345	406	II
X9376538140_R04C01	3212	550	II

3. Fill the following table: what are the Red and Green fluorescence for the 41620492 address? Optional: check from the manifest file if the address corresponds to a Type I or a Type II probe and, in case of Type I probe, report its color.

```
#head(Red)
red_address <- Red[rownames(Red) == "41620492",]
red_address
#head(Green)
green_address <- Green[rownames(Green) == "41620492",]
green_address


load('/home/maria/Desktop/DRD2/Illumina450Manifest.RData')
type<-as.character(Illumina450Manifest[Illumina450Manifest$
    AddressA_ID == "41620492",]$Infinium_Design_Type)
type
allcolors <- as.data.frame(cbind(t(red_address),t(green_address)))
allcolors</pre>
```

4. Create the object MSet.raw .

- 5. Perform the following quality checks and provide a brief comment to each step:
  - QCplot:

```
qc <- getQC(MSet.raw)
plotQC(qc)</pre>
```


getQC minfi: Estimates sample-specific quality control (QC) for methylation data. It gets the value of a DataFrame with two columns: mMed and uMed which are the chipwide medians of the Meth and Unmeth channels. Then it provides a plot where it is possible to observe that good samples cluster together, while failed samples tend to separate and have lower median intensities.



• check the intensity of negative controls using minfi:

Internal control probes can be used to assess the quality of different sample preparation steps (bisulfite conversion, hybridization, etc.). The values of these control probes are stored in the initial RGChannelSet and can be plotted by using the function control-StripPlot and by specifying the control probe type that in this case are Negative controls

### Control: NEGATIVE



• Calculate detection pValues; for each sample, how many probes have a detection p-value higher than the threshold assigned to each student?

```
### Detection p-value
##########################

?detectionP

detP <- detectionP(RGset)

save(detP,file="detP.RData")

str(detP)

dim(detP)

head(detP)

failed <- detP>0.05
head(failed)
dim(failed)
table(failed)

summary(failed) # In this way we can now how many probes
failed for each sample
```

we used the detectionP function of minfi: this function detect the p-values, for all probed genomic positions in every sample. Small p-values indicate a good position. Positions

SAMPLE	Failed positions
9344737127_R01C02	499
9344737127_R03C02	428
9344737127_R06C01	480
9344737127_R02C02	463
9376538140_R02C01	441
9376538140_R05C01	420
9376538140_R01C01	523
9376538140_R04C01	97

with non-significant p-values (here > 0.05) should not be trusted. So we identified failed positions defined as both the methylated and unmethylated channel reporting background signal levels.

• Create a vector of bad probes to be removed because they have a detection pValue higher that the assigned threshold (see above) in more than 1% of the samples.

```
# We can remove the probes having > 1 % of samples with a
detection p-value greater than 0.05

means_rows <- rowMeans(failed)
head(means_rows)
length(means_rows)

probes_to_be_removed <- means_rows[means_rows>0.01]
length(probes_to_be_removed)
head(probes_to_be_removed)
names_probes_to_be_removed <- names(probes_to_be_removed)
head(names_probes_to_be_removed)
str(names_probes_to_be_removed)
save(names_probes_to_be_removed, file="names_probes_to_be_
removed.RData")
```

6. Calculate raw beta and M values and plot the densities of mean methylation values, dividing the samples in DS and WT .

Beta value is the percentage of methylation at the given CpG site; it is a continuous variable between 0 (absent methylation) and 1(completely methylated)

$$\beta = \frac{M}{M+U} \tag{1}$$

The M value is a continuous variable which can in principle take on any value on the real line.

$$M = \log_2 \frac{M}{U} \tag{2}$$


```
#subset the beta and M values matrixes in order to retain DS or
    WT subjects
DS <- targets[targets$Group=="DS",]
DS <- droplevels(DS)
RGset_DS <- read.metharray.exp(targets = DS)
MSet_DS.raw <- preprocessRaw(RGset_DS)

WT <- targets[targets$Group=="WT",]
WT <- droplevels(WT)
RGset_WT <- read.metharray.exp(targets = WT)
MSet_WT.raw <- preprocessRaw(RGset_WT)</pre>
```

• Beta value:

```
#Density plots of Beta:
2 # I want to plot the density distribution of mean beta for each
3 # First, I have to calculate the mean of each row of the
     dataframe
4 # apply <- works on a dataframe or a matrix and returns a vector
5 ?apply
6 beta_DS <- getBeta(MSet_DS.raw)
p beta_WT <- getBeta(MSet_WT.raw)</pre>
9 mean_beta_DS <- apply(beta_DS,1,mean) # 1 indicates that I apply
      the functiont to each row
mean_beta_WT <- apply(beta_WT,1,mean)</pre>
d_mean_beta_DS <- density(mean_beta_DS, na.rm=T)
d_mean_beta_WT <- density(mean_beta_WT,na.rm=T)
par (mfrow=c(1,2))
plot(d_mean_beta_WT,col=2,main="Density of Mean Beta-values",
     xlab="Beta")
15 lines(d_mean_beta_DS,col=7)
16 legend('topright',pch=c(16,16),col=c(2,7),legend=c("WT","DS"))
```

## **Density of Mean Beta-values**




#### • M value:

```
#Density plots of M values :
# I want to plot the density distribution of mean M values for
each CpG

Mvalue_DS <- getM(MSet_DS.raw)
Mvalue_WT <- getM(MSet_WT.raw)
mean_M_DS <- apply(Mvalue_DS,1,mean)
mean_M_WT <- apply(Mvalue_WT,1,mean)
d_mean_M_DS <- density(mean_M_DS,na.rm=T)
d_mean_M_WT <- density(mean_M_WT,na.rm=T)
par(mfrow=c(1,2))
plot(d_mean_M_WT,col=3,main="Density of Mean M-values",xlab="
Beta")
lines(d_mean_M_DS,col=6)
legend('topright',pch=c(16,16),col=c(3,6),legend=c("WT","DS"))</pre>
```


# **Density of Mean M-values**



7. Compare raw data with normalized data using preprocessQuantile function. Produce a plot with 6 panels (2 rows, 3 columns) in which, for both raw and normalized data, you show the density plots of beta mean values according to the chemistry of the probes, the density plot of beta standard deviation values according to the chemistry of the probes and the boxplot of beta values. Provide a short comment regarding the changes you observe.

```
| load('/home/maria/Desktop/DRD2/Illumina450Manifest_clean.RData')
3 dfI <- Illumina450Manifest_clean[Illumina450Manifest_clean$
     Infinium_Design_Type == "I",]
4 dfI <- droplevels(dfI)
5 dim (dfI)
6 str(dfI)
 dfII <- Illumina450Manifest_clean[Illumina450Manifest_clean$
     Infinium_Design_Type=="II",]
 dfII <- droplevels(dfII)</pre>
10
 # I subset the dataframe beta in order to retain only the rows
     whose name is in the first column of dfI
12 beta_I <- beta[rownames(beta) %in% dfI$IlmnID,]</pre>
13 dim(beta_I)
 # I subset the dataframe beta in order to retain only the rows
     whose name is in the first column of dfII
beta_II <- beta[rownames(beta) %in% dfII$IlmnID,]
16 dim(beta_II)
17
18 # I want to compare the density distributions of type I and
     type II probes
 # first, I have to calculate the mean of beta for each row
20 mean_of_beta_I <- apply(beta_I,1,mean)</pre>
mean_of_beta_II <- apply(beta_II,1,mean)</pre>
23 d_mean_of_beta_I <- density(mean_of_beta_I,na.rm=T)
d_mean_of_beta_II <- density(mean_of_beta_II,na.rm=T)
25
 plot(d_mean_of_beta_I,main="Mean Of Raw Beta-values",col=2,xlab=
     "mean of beta")
27 lines(d_mean_of_beta_II,col=4)
28 legend('topright', pch=c(16,16), col=c(2,4), legend=c("beta_I","
     beta_II"))
29 #####
30 # For Raw data we have just prepared the density plot of mean
     betas...let's prepare the density plot of standard deviations
 sd_of_beta_I <- apply(beta_I,1,sd)</pre>
sd_of_beta_II <- apply(beta_II,1,sd)</pre>
d_sd_of_beta_I <- density(sd_of_beta_I,na.rm=T)
d_sd_of_beta_II <- density(sd_of_beta_II,na.rm=T)</pre>
36 plot(d_sd_of_beta_I,main="Standard deviation of Raw Beta-values"
     , col=2, xlab="beta standard deviation")
```

```
37 lines(d_sd_of_beta_II,col=4)
38 legend('topright', pch=c(16,16), col=c(2,4), legend=c("beta_I","
     beta_II"))
39
 boxplot(beta, main="Raw Beta-values", las=2, cex.axis=0.5)
40
41
 42
43
 # preprocessQuantile (minfi; input: RGset or Mset; output:
44
     GenomicRatioSet)
?preprocessQuantile
46 RGset
47 preprocessQuantile_results <- preprocessQuantile(RGset)
48 str(preprocessQuantile_results)
49 class (preprocessQuantile_results)
50 preprocessQuantile_results
51 beta_preprocessQuantile <- getBeta(preprocessQuantile_results)</pre>
52 head(beta_preprocessQuantile)
53
54 beta_preprocessQuantile_I <- beta_preprocessQuantile[rownames(</pre>
     beta_preprocessQuantile) %in% dfI$IlmnID,]
 beta_preprocessQuantile_II <- beta_preprocessQuantile[rownames(</pre>
     beta_preprocessQuantile) %in% dfII$IlmnID,]
56 mean_of_beta_preprocessQuantile_I <- apply(beta_
     preprocessQuantile_I,1,mean)
mean_of_beta_preprocessQuantile_II <- apply(beta_
     preprocessQuantile_II,1,mean)
58 d_mean_of_beta_preprocessQuantile_I <- density(mean_of_beta_
     preprocessQuantile_I, na.rm=T)
59 d_mean_of_beta_preprocessQuantile_II <- density(mean_of_beta_
     preprocessQuantile_II, na.rm=T)
60 sd_of_beta_preprocessQuantile_I <- apply(beta_preprocessQuantile
     _I,1,sd)
61 sd_of_beta_preprocessQuantile_II <- apply(beta_
     preprocessQuantile_II,1,sd)
62 d_sd_of_beta_preprocessQuantile_I <- density(sd_of_beta_
     preprocessQuantile_I, na.rm=T)
63 d_sd_of_beta_preprocessQuantile_II <- density(sd_of_beta_
     preprocessQuantile_II, na.rm=T)
_{64} par (mfrow=c(1,3))
65 jpeg('norm1.jpg')
66 plot(d_mean_of_beta_preprocessQuantile_I,na.rm=T,main="Mean of
     normalized beta values", col=2, xlab="mean of beta")
67 lines(d_mean_of_beta_preprocessQuantile_II,col=4)
 legend('topright', pch=c(16,16), col=c(2,4), legend=c("beta_I","
     beta_II"))
```



From these plots it can be appreciated how the normalization procedure removed technical variation between different array and different samples. As we can see from the plot, after applying preprocessQuantile normalization, the distributions of -values appear to look more similar than Raw Beta-values.


8. Consider the normalized beta and M values and remove the probes previously defined as bad according to the detection pValue.


```
beta_preprocessQuantile <- getBeta(preprocessQuantile_results)
M_preprocessQuantile <- getM(preprocessQuantile_results)


summary(beta_preprocessQuantile)
summary(M_preprocessQuantile)
beta_filt_preprocessQuantile <- beta_preprocessQuantile[!
rownames(beta_preprocessQuantile) %in% names_probes_to_be_</pre>
```


9. Check beta and M matrixes generated in step 8 for homo/heteroschedasticity; comment the plot. Optional: plot the lowess line.

```
# HOMOSCEDASTICITY and HETEROSCEDASTICITY
 beta_filt_preprocessQuantile_mean <- apply(beta_filt_</pre>
     preprocessQuantile,1, mean)
 beta_filt_preprocessQuantile_stdev <- apply(beta_filt_</pre>
     preprocessQuantile,1, sd)
6 M_filt_preprocessQuantile_mean <- apply(M_filt_
     preprocessQuantile,1, mean)
7 M_filt_preprocessQuantile_stdev <- apply(M_filt_
     preprocessQuantile,1, sd)
 jpeg('betameandsd3.jpg')
 plot(beta_filt_preprocessQuantile_mean, beta_filt_
     preprocessQuantile_stdev)
10 dev.off()
jpeg('Mmeandsd3.jpg')
plot(M_filt_preprocessQuantile_mean, M_filt_preprocessQuantile_
     stdev)
dev.off()
15 par (mfrow=c(1,2))
jpeg('LOWESS1.jpg')
17 smoothScatter(beta_filt_preprocessQuantile_mean, beta_filt_
     preprocessQuantile_stdev)
18 lines(lowess(beta_filt_preprocessQuantile_mean, beta_filt_
     preprocessQuantile_stdev), col="red") # lowess carries out a
     locally weighted regression of y on x
dev.off()
jpeg('LOWESS2.jpg')
smoothScatter(M_filt_preprocessQuantile_mean, M_filt_
     preprocessQuantile_stdev)
22 lines(lowess(M_filt_preprocessQuantile_mean, M_filt_
     preprocessQuantile_stdev), col="red") # lowess line (x,y)
23 dev.off()
```










This plots show Beta-value with a significant heteroscedasticity and M-values are approximately more homoscedastic.

10. Perform a PCA on the beta matrix generated in step 8 to check for batch effects. Comment the plot.

### Normalized beta values PCA



The distances among samples in a pca scores plot are related to their similarity with respect to the pattern described by the model and so this type of plots can be used for the identification of outliers and batch effects. Batch effects are technical sources of variation that have been added to the samples during handling. Here, from the obtained plot no obvious batch effect is observed.

11. Using the beta matrix generated in step 8, identify differentially methylated probes between group A and group B using Anova correcting for Female .

```
#####11####

pheno <- read.csv("/home/maria/Desktop/DRD2/Input_data/
SampleSheet.csv",sep=",")
```

#### Results:

> head(final)

```
X9344737127_R01C02 X9344737127_R03C02 X9344737127_R06C01 X9344737127_R02C02
cg13869341
               0.73914779
                                0.76950459
                                                     0.6995062
                                                                     0.7184184
cg14008030
                 0.69185593
                                 0.64592974
                                                     0.6885158
                                                                      0.6031361
                                 0.11699326
cg12045430
                0.13474718
                                                   0.1004484
                                                                      0.1218535
                                 0.28879485
cg20826792
                 0.17078769
                                                    0.2947075
                                                                      0.2819566
                                 0.09592991
cg00381604
                 0.05846309
                                                    0.1045477
                                                                      0.1000693
                                 0.51681548
cg20253340
                 0.42958351
                                                     0.4396129
                                                                      0.4760213
         X9376538140_R02C01 X9376538140_R05C01 X9376538140_R01C01 X9376538140_R04C01
cg13869341
               0.78299289 0.77009158
                                                   0.7751171
                                                                      0.8177119
cg14008030
                0.63207688
                                  0.56279678
                                                    0.6861439
                                                                      0.6460894
cg12045430
                0.14349114
                                 0.12472814
                                                   0.4254814
                                                                     0.1188433
                 0.22464011
                                 0.20651747
                                                   0.3439327
cg20826792
                                                                     0.2847663
                                0.08688974
cg00381604
                 0.09964109
                                                   0.3128018
                                                                      0.1086204
cg20253340
                 0.46439970
                                0.37458736
                                                    0.4502786
                                                                      0.6004477
         pValuesAnova
cg13869341 0.68847287
cg14008030 0.40103860
cg12045430 0.46427208
cg20826792 0.05315406
cg00381604 0.24584491
cg20253340 0.38870914
```

12. Apply multiple test correction and set a significant threshold of 0.05. How many probes do you identify as differentially methylated considering nominal pValues? How many after Bonferroni correction? How many after BH correction?


```
###12###
length(pValuesAnova[pValuesAnova<=0.05])
pValuesAnovaBonf <- p.adjust(pValuesAnova, "bonferroni")
length(pValuesAnovaBonf[pValuesAnovaBonf<=0.05])
pValuesAnovaBH <- p.adjust(pValuesAnova, "BH")</pre>
```


```
6 length (pValuesAnovaBH [pValuesAnovaBH <= 0.05])
```

Probes resulted to be differentially methylated were 36693. But after both Bonferroni correction and BH correction they resulted to be 0.

13. Produce a Manhattan plot and a volcano plot of your data.

```
# VOLCANO PLOTS
  ****************
 #For each probe of the microarray, we have to calculate the
     delta between mean group A values and mean group B values
 beta <- final[,1:8]
6 beta_groupA <- beta[,pheno$Group=="A"]
 mean_beta_groupA <- apply(beta_groupA,1,mean)</pre>
 beta_groupB <- beta[,pheno$Group=="B"]</pre>
 mean_beta_groupB <- apply(beta_groupB,1,mean)</pre>
10 delta <- mean_beta_groupB-mean_beta_groupA
11 head (delta)
12
13 toVolcPlot <- data.frame(delta, -log10(final$pValuesAnova))
14 head (to VolcPlot)
plot(toVolcPlot[,1], toVolcPlot[,2])
16
 # MANHATTAN PLOTS
 18
  #manhattan plot shows significant methylation sites in the
     genome
20 install.packages("gap")
21 library (gap)
22 # First we have to annotate our dataframe, that is add genome
     annotation information for each cpg probe
23 final <- data.frame(rownames(final), final)
24 colnames(final)[1] <- "IlmnID"
 annotated <- merge(final, Illumina450Manifest_clean,by="IlmnID")
26 save (annotated, file="final_annotated.RData")
 db <- data.frame(annotated$CHR, annotated$MAPINFO, annotated$
     pValuesAnova)
29 levels (db$annotated.CHR)
 db$annotated.CHR <- factor(db$annotated.CHR,levels=c("1","2","3"
     ,"4","5","6","7","8","9","10","11","12","13","14","15","16","
     17","18","19","20","21","22","X","Y"))
31 levels (db$annotated.CHR)
32 palette <- rainbow(24)
mhtplot(db,control=mht.control(colors=palette))
```





 $14.\ \,$  Produce an heatmap of the top 100 differentially mehtylated probes.

```
# HEATMAP
```

```
install.packages("gplots")
library(gplots)

# Heatmap on several probes is computationally demonading. We use just the top 100 most significant CpG probes
annotated <-annotated[order(annotated$pValuesAnova),]

matrix=as.matrix(annotated[1:100,2:9])
colorbar <- c("green", "green", "orange", "orange", "green", "orange", "orange", "orange")
heatmap.2(matrix,col=terrain.colors(100),Rowv=T,Colv=T, dendrogram="both",key=T,ColSideColors=colorbar,density.info="none",trace="none",scale="none",symm=F)</pre>
```

